

Recombinant Sheep Interleukin-8(IL8)

Catalog Number: CSB-EP011671SH



Product Name:	Recombinant Sheep Interleukin-8(IL8)
Alternative names:	C-X-C motif chemokine 8
Catalog Number:	CSB-EP011671SH
Relevance :	IL-8 is a chemotactic factor that attracts neutrophils, basophils, and T-cells, but not monocytes. It is also involved in neutrophil activation. It is released from several cell types in response to an inflammatory stimulus.
Mol. Weight:	14kd
Product Info :	His-tagged
Source:	E.coli derived
Image:	
Purity:	>90%(SDS-PAGE)
Storage Buffer:	20mM Tris-HCl, 0.5M NaCl, pH 8.0,50% glycerol
Storage :	Store at -20℃, for extended storage, conserve at -20℃ or -80℃.
Notes :	Repeated freezing and thawing is not recommended. Store working aliquots at 4℃ for up to one week.
AA sequence:	AVLSRMSTELRCQCICKTHSTPFHPKFIKELRVIESGPHCENSEIIVKLTNGKEVCLDPKE KVVQKVVAFLKRAEKQDP
References:	"Sequencing of the ovine interleukin-8-encoding cDNA using the polymerase chain reaction." Legastelois I., Greenland T., Arnaud P., Mornex J.F., Cordier G. Gene 150:367-369(1994)