




Recombinant E.coli 30S ribosomal protein S2

Catalog Number: CSB-RP086444Ba

Product Name:	Recombinant E.coli 30S ribosomal protein S2
Catalog Number:	CSB-RP086444Ba
Relevance :	Part of the 30S ribosomal subunit. Some nascent polypeptide chains are able to cross-link to this protein in situ.
Mol. Weight:	53kD
Product Info :	GST-tagged
Source:	E.coli derived
Images	
Purity:	>90%(SDS-PAGE)
Storage Buffer:	PBS, pH 7.4, 50% glycerol
Storage :	Store at -20℃, for extended storage, conserve at -20℃ or -80℃.
Notes :	Repeated freezing and thawing is not recommended. Store working aliquots at 4℃ for up to one week.
AA sequence:	ATVSMRDMLKAGVHFGHQTRYWNPCKMPFIFGARNKVHIINLEKTVPMFNEALAE LNKIASRK GKILFVGTKRAASEAVKDAALSCDQFFVNHRWLGGMLTNWKTVRQSIKRLKDLETQSQDGTG DKLTKEALMRTRELEKLENSLGGIKDMGGLPDALFVIDADHEHIAIKEANNLGIPVFAIVDTNSD PDGVDFVIPGNDDAIRAVTLYLGAVAATVREGRSQDLASQAEESFVEAE
References:	"Organization and nucleotide sequence of a new ribosomal operon in Escherichia coli containing the genes for ribosomal protein S2 and elongation factor Ts." An G., Bendiak D.S., Mamelak L.A., Friesen J.D. Nucleic Acids Res. 9:4163-4172(1981)