




Recombinant Mouse Vascular endothelial growth factor C (Vegfc)

Product Code	CSB-YP025835MO
Relevance	Growth factor active in angiogenesis, and endothelial cell growth, stimulating their proliferation and migration and also has effects on the permeability of blood vessels. May function in angiogenesis of the venous and lymphatic vascular systems during embryogenesis, and also in the maintenance of differentiated lymphatic endothelium in adults. Binds and activates VEGFR-2 (KDR/FLK1) and VEGFR-3 (FLT4) receptors.
Abbreviation	Recombinant Mouse Vegfc protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P97953
Alias	Flt4 ligand ;Flt4-L Vascular endothelial growth factor-related protein ;VRP
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	AHYNTEILKSIDNEWWRKTQCMPREVCIDVGKEFGAATNTFFKPPCVSVYRCGG CCNSEGLQCMNTSTGYLSKTLFEITVPLSQGPKPVTISFANHTSCRCMSKLDV YRQVHSIIRR
Research Area	Others
Source	Yeast
Target Names	Vegfc
Protein Names	Recommended name: Vascular endothelial growth factor C Short name= VEGF-C Alternative name(s): Flt4 ligand Short name= Flt4-L Vascular endothelial growth factor-related protein Short name= VRP
Expression Region	108-223aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	15.0kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.