



# Recombinant Mouse Sex-determining region Y protein (Sry), partial

<b>Product Code</b>	CSB-YP022699MO
<b>Relevance</b>	Transcriptional regulator that controls a genetic switch in male development. It is necessary and sufficient for initiating male sex determination by directing the development of supporting cell precursors (pre-Sertoli cells) as Sertoli rather than granulosa cells. In male adult brain involved in the maintenance of motor functions of dopaminergic neurons . Involved in different aspects of gene regulation including promoter activation or repression. SRY HMG box recognizes DNA by partial intercalation in the minor groove. Promotes DNA bending. Also involved in pre-mRNA splicing . Binds to the DNA consensus sequence 5'-[AT]AACAA[AT]-3'.1 Publication
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q05738
<b>Alias</b>	Testis-determining factor
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	MEGHVKRPMNAFMVWSRGERHKLAAQQNPSMQNTEISKQLGCRWKSLEAE KRPFFQEAQRLKILHREKYPNYKYQPHHRAKVSQSRSGILQPAVASTKLYNLLQ WDRNPHAITYRQDWSRAAHLYSKNQQSFYWQPVDIPTGHLQ
<b>Research Area</b>	Others
<b>Source</b>	Yeast
<b>Target Names</b>	Sry
<b>Protein Names</b>	Recommended name: Sex-determining region Y protein Alternative name(s): Testis-determining factor
<b>Expression Region</b>	1-144aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.