




Recombinant Human Alcohol dehydrogenase class-3 (ADH5)

Product Code	CSB-YP001357HU
Relevance	Class-III ADH is rarkably ineffective in oxidizing ethanol, but it readily catalyzes the oxidation of long-chain primary alcohols and the oxidation of S-(hydroxymethyl) glutathione.
Abbreviation	Recombinant Human ADH5 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P11766
Alias	Alcohol dehydrogenase 5Alcohol dehydrogenase class chi chain;Alcohol dehydrogenase class-IIIGlutathione-dependent formaldehyde dehydrogenase (EC:1.1.1.-) ;FALDH ;FDH ;GSH-FDHS-(hydroxymethyl)glutathione dehydrogenase (EC:1.1.1.284)
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	ANEVIKCKAAVAWEAGKPLSIEEIEVAPPKAHEVRIKIIATAVCHTDAYTLSGADP EGCFPVILGHEGAGIVESVGEGVTKLKAGDTVIPLYIPQCGECKFCLNPKTNLC QKIRVTQGKGLMPDGTSRFTCKGKTIHYMGSTSTFSEYTVVADISVAKIDPLAP LDKVCLLGCGISTGYGAAVNTAKLEPGSVCAVFGGLGGVGLAVIMGCKVAGASR IIGVDINKDKFARAKEFGATECINPQDFSKPIQEVLIEMTDGGVDYSFECIGNVK VMRAALEACHKGWGVSVVVGVAASGEEIATRPFLVTGRTWKGTAFGGWKS VESVPKLVSEYMSKKIKVDEFVTHNLSFDEINKAFELMHSGKSIKRTVVKI
Research Area	Metabolism
Source	Yeast
Target Names	ADH5
Protein Names	Recommended name: Alcohol dehydrogenase class-3 EC= 1.1.1.1 Alternative name(s): Alcohol dehydrogenase 5 Alcohol dehydrogenase class chi chain Alcohol dehydrogenase class-III Glutathione-dependent formaldehyde dehydrogenase
Expression Region	2-374aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	41.6kDa

**Protein Length**

Full Length of Mature Protein

Image

(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.