



# Recombinant Sheep Interferon gamma (IFNG)

<b>Product Code</b>	CSB-YP011050SH
<b>Relevance</b>	Produced by lymphocytes activated by specific antigens or mitogens. IFN-gamma, in addition to having antiviral activity, has important immunoregulatory functions. It is a potent activator of macrophages, it has antiproliferative effects on transformed cells and it can potentiate the antiviral and antitumor effects of the type I interferons.
<b>Abbreviation</b>	Recombinant Sheep IFNG protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P17773
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Ovis aries (Sheep)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	QGPFFKEIENLKEYFNASNPDVAKGGPLFSEILKNWKEESDKKIIQSQIVSFYFK LFENLKDQVIQRSMDIKQDMFQKFLNGSSEKLEDFKRLIQIPVDDLQIQRKAI NELIKVMNDLSPKSNLRKRKRSQNLFRGRRASM
<b>Research Area</b>	Others
<b>Source</b>	Yeast
<b>Target Names</b>	IFNG
<b>Protein Names</b>	Recommended name: Interferon gamma Short name= IFN-gamma
<b>Expression Region</b>	24-166a
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	18.9kDa
<b>Protein Length</b>	Full Length of Mature Protein
<b>Image</b>	

