




Recombinant Mouse Glutathione S-transferase P 1 (Gstp1)

| | |
|--------------------------|---|
| Product Code | CSB-YP009989MO |
| Relevance | Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles. Can metabolize 1-chloro-2,4-dinitrobenzene. Regulates negatively CDK5 activity via p25/p35 translocation to prevent neurodegeneration . |
| Abbreviation | Recombinant Mouse Gstp1 protein |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P19157 |
| Alias | GST YF-YFGST class-piGST-piBPreadipocyte growth factor |
| Product Type | Recombinant Protein |
| Immunogen Species | Mus musculus (Mouse) |
| Purity | ≥ 90% as determined by SDS-PAGE. |
| Sequence | PPYTIVYFPVVRGRCEAMRMLLADQGGQSWKEEVVTIDTWMQGLLKPTCLYGQL PKFEDGDLTLYQSNAILRHLGRSLGLYGKNQREAAQMDMVNDGVEDLRGKYV TLIYTNYENGKNDYVKALPGHLKPFETLLSQNQGGKAFIVGDQISFADYNLLDL LLIHQVLPAGCLDNFPLLSAYVARLSARPKIKAFLLSSPEHVNRPINGNGKQ |
| Research Area | Others |
| Source | Yeast |
| Target Names | Gstp1 |
| Protein Names | Recommended name: Glutathione S-transferase P 1 Short name= Gst P1 EC= 2.5.1.18 Alternative name(s): GST YF-YF GST class-pi GST-piB Preadipocyte growth factor |
| Expression Region | 2-210aa |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | N-terminal 6xHis-tagged |
| Mol. Weight | 25.5kDa |
| Protein Length | Full Length of Mature Protein |
| Image | |



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.