



# Recombinant Rat Perforin-1 (PRF1), partial

<b>Product Code</b>	CSB-MP018668RA
<b>Relevance</b>	Plays a key role in secretory granule-dependent cell death, and in defense against virus-infected or neoplastic cells. Can insert into the mbrane of target cells in its calcium-bound form, oligomerize and form large pores. Promotes cytolysis and apoptosis of target cells by facilitating the uptake of cytotoxic granzymes.
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P35763
<b>Alias</b>	Cytolysin;Lymphocyte pore-forming protein
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	<p>           ATRSECKQNHKFVPGVWAAGEGVDVTTLRSSSFPVNTGKFLRPDRTCTLCK            NALMNDGIQRLPVAIAHWRPHGSHCQRNVATTKVSSTEGVAREAAANINNDW            RAGLDVNPKEANVHVSVAGSHSKIANFAAEKAHQDQYNFNTDTVECRMYSF            RLAQKPPLHPDFRKALKNLPHNFNSSTEHAYRRLISSYGTHFITAVDLGGRVSV            LTALRTCQLTLDGLTADEVGDCLSVEAQVSIGAQASVSSEYKACEEKKKQHKIA            TSFHQTYRERHVEVLGGPLDSSNDLLFGNQATPEHFSTWIASLPTRPDVVDYS            LEPLHILLEDSDPKREALRQAISHYVMSRARWRDCNRPCRAGQHKSSRDSCQ            CVCQDSNVTNQDCCPRQRGLAKLMVRNFQAKGLWGDYITSTDAYLKVFFGG            QEIRTGVVWNNNHPSWSDKMDFGNVLLSTGGPLRVQVWDADNGWDDDLLG            TCDKSPKSGFHEVNCPLNHGSIKFIYQANCLPDLTGETCLEYAPQGLLGDPRG            NRSGAVW         </p>
<b>Research Area</b>	Others
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Prf1
<b>Protein Names</b>	Recommended name: Perforin-1 Short name= P1 Alternative name(s): Cytolysin Lymphocyte pore-forming protein
<b>Expression Region</b>	25-554aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.