




Recombinant Mouse Superoxide dismutase [Cu-Zn] (Sod1)

Product Code	CSB-YP022397MO
Relevance	Destroys radicals which are normally produced within the cells and which are toxic to biological systems.
Abbreviation	Recombinant Mouse Sod1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P08228
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	AMKAVCVLKGDPVQGTHFEQKASGEPVVLSGQITGLTEGQHGFHVHQYGD NTQGCTSAGPHFNPHSKKHGGPADEERHVGDLGNVTAGKDGVANVSIEDRVI SLSGEHSIIGRTMVVHEKQDDLKGGNEESTKTGNAGSRLACGVIGIAQ
Research Area	Others
Source	Yeast
Target Names	Sod1
Protein Names	Recommended name: Superoxide dismutase [Cu-Zn] EC= 1.15.1.1
Expression Region	2-154aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	17.8kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.