




Recombinant Human Corticosteroid-binding globulin (SERPINA6)

Product Code	CSB-YP021062HU
Relevance	Major transport protein for glucocorticoids and progestins in the blood of almost all vertebrate species.
Abbreviation	Recombinant Human SERPINA6 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P08185
Alias	Serpin A6Transcortin
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MDPNAAYVNMSNHHRGLASANVDFAFSLYKHLVALSPKKNIFISPVISMALAM LSLGTGCGHTRAQLLQGLGFNLTERSETEIHQGFQHLHQLFAKSDTSLEMTMGN ALFLDGSLELLESFSADIKHYESEVLAMNFQDWATASRQINSYVKNKTQGGKIV DLFSGLDSPAILVLVNYIFFKGTWTQPFDLASTREENFYVDETTVVKVPMMLQS STISYLHDSLEPCQLVQMNYVGNGTVFFILPDKGKMNTVIAALSRDTINRWSAG LTSSQVDLYIPKVTISGVYDLGDVLEEMGIADLFTNQANFSRITQDAQLKSSKVV HKAVLQLNEEGVDTAGSTGVTLNLTSPKPIILRFNQPFIMIFDHFTWSSLFLARV MNPV
Research Area	Transport
Source	Yeast
Target Names	SERPINA6
Protein Names	Recommended name: Corticosteroid-binding globulin Short name= CBG Alternative name(s): Serpin A6 Transcortin
Expression Region	23-405aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	44.6kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.