




# Recombinant Human Interferon alpha-10 (IFNA10)

<b>Product Code</b>	CSB-YP011033HU
<b>Relevance</b>	Produced by macrophages, IFN-alpha have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase.
<b>Abbreviation</b>	Recombinant Human IFNA10 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P01566
<b>Alias</b>	Interferon alpha-6L Interferon alpha-C ;LeIF C
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	CDLPQTHSLGNRRALILLGQMGRISPFSLCKDRHDFRIPQEEFDGNQFQKAQA ISVLHEMIQQTFNLFSTEDSSAAWEQSLLEKFSTELYQQLNDLEACVIQEVGVE ETPLMNEDSILAVRKYFQRITLYLIERKYSPLCAWEVVRAEIMRSLSFSTNLQKRL RRKD
<b>Research Area</b>	Immunology
<b>Source</b>	Yeast
<b>Target Names</b>	IFNA10
<b>Protein Names</b>	Recommended name: Interferon alpha-10 Short name= IFN-alpha-10 Alternative name(s): Interferon alpha-6L Interferon alpha-C Short name= LeIF C
<b>Expression Region</b>	24-189aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	21.4kDa
<b>Protein Length</b>	Full Length of Mature Protein
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.