



# Recombinant Human Fibroblast growth factor 9 (FGF9)

<b>Product Code</b>	CSB-BP008638HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P31371
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	LGEVGN Y FGVQDAVPFG NVPVLPVDSP VLLSDHLGQS EAGGLPRGPA VTDL D H L K G I L R R R Q L Y C R T G F H L E I F P N G T I Q G T R K D H S R F G I L E F I S I A V G L V S I R G V D S G L Y L G M N E K G E L Y G S E K L T Q E C V F R E Q F E E N W Y N T Y S S N L Y K H V D T G R R Y Y V A L N K D G T P R E G T R T K R H Q K F T H F L P R P V D P D K V P E L Y K D I L S Q S
<b>Source</b>	Baculovirus
<b>Target Names</b>	FGF9
<b>Protein Names</b>	Recommended name: Fibroblast growth factor 9 Short name= FGF-9 Alternative name(s): Glia-activating factor Short name= GAF Heparin-binding growth factor 9 Short name= HBGF-9
<b>Expression Region</b>	4-208
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This protein is a member of the fibroblast growth factor (FGF) family. FGF family members possess broad mitogenic and cell survival activities, and are involved in a variety of biological processes, including embryonic development, cell growth, morphogenesis, tissue repair, tumor growth and invasion. This protein was isolated as a secreted factor that exhibits a growth-stimulating effect on cultured glial cells. In nervous system, this protein is produced mainly by neurons and may be important for glial cell development. Expression of the mouse homolog of this gene was found to be dependent on Sonic hedgehog (Shh) signaling. Mice lacking the homolog gene displayed a male-to-female sex reversal phenotype, which suggested a role in testicular embryogenesis.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.