



Recombinant Human RAC-beta serine/threonine-protein kinase (AKT2), partial

Product Code	CSB-BP001555HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P31751
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	NEVSVIKEG WLHKRGEYIK TWRPRYFLLK SDGSFIGYKE RPEAPDQTLPLNNSVAEC QLMKTERPRP NTFVIRCLQW TTVIERTFHV DSPDEREEWMRAIQMVANSL KQRAPGEDPM DYKCGSPSDS STTEEMEVAV SKARAKVTMNDFDYLKLLGK GTFGKVLVR EKATGRYYAM KILRKEVIA KDEVAHTVTE SRVLQNTRHP FLTALKYAFQ THDRLCFVME YANGGELFFH LSRERVFTEERARFYGAEIV SALEYLHSRD VVYRDIKLEN LMLDKDGHK ITDFGLCKEGISDGATMKTFCGTPEYLAPE VLEDNDYGRA VDWWGLGVVMYEMMCGRLPF YNQDHERLFE LILMEEIRFP RTLSPEAKSL LAGLLKKDPKQRLGGGPSDA KEVMEHRFFL SINWQDVVQK KLLPPFKPQV TSEVDTRYFDDEFTAQSITI TPPDRYDSLGLLELDQRTHF PQFSYSASIR E
Source	Baculovirus
Target Names	AKT2
Protein Names	Recommended name: RAC-beta serine/threonine-protein kinase EC=2.7.11.1Alternative name(s): Protein kinase Akt-2 Protein kinase B beta Short name= PKB beta RAC-PK-beta
Expression Region	2-481
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Target Details	This gene is a putative oncogene encoding a protein belonging to a subfamily of serine/threonine kinases containing SH2-like (Src homology 2-like) domains. The gene was shown to be amplified and overexpressed in 2 of 8 ovarian carcinoma cell lines and 2 of 15 primary ovarian tumors. Overexpression contributes to the malignant phenotype of a subset of human ductal pancreatic cancers. The encoded protein is a general protein kinase capable of phosphorylating several known proteins.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.