



# Recombinant Pig Bone sialoprotein 2 (IBSP)

<b>Product Code</b>	CSB-MP010945PI
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P31936
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Sus scrofa (Pig)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	FSMKNFHRRRA KLEDPEENGV FKYRPRYYLY KHAYFYPLK RFPVQSSSDS SEENGN GDSS EEEEEEEENS NEEENNEENE DSDGNEDEDS EAENITLSTT TLGYGGDVTP GTASIGLAAL QLPKKAGDIG KKSAKEEESD EEEEEEEENE ENEA EVDNE QGTNGTSTNS TEVD SGN GHS GGDNGEEGDQ ESVTEAQGTT VAGEQDNGGA KTTTSPNGGL EPTPPPQDIS GTTLPPSGKT TTPEYEGEYE QTGAHEYDNG YEIYSENGE PRGDSYRAYE DEYSYYKGRS YNSYGGHDYY
<b>Source</b>	Mammalian cell
<b>Target Names</b>	IBSP
<b>Protein Names</b>	Recommended name: Bone sialoprotein 2 Alternative name(s): Bone sialoprotein II Short name= BSP II Cell-binding sialoprotein Integrin-binding sialoprotein
<b>Expression Region</b>	1-300
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	This protein is a major structural protein of the bone matrix. It constitutes approximately 12% of the noncollagenous proteins in human bone and is synthesized by skeletal-associated cell types, including hypertrophic chondrocytes, osteoblasts, osteocytes, and osteoclasts. The only extraskelatal site of its synthesis is the trophoblast. This protein binds to calcium and hydroxyapatite via its acidic amino acid clusters, and mediates cell attachment through an RGD sequence that recognizes the vitronectin receptor.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.