



# Recombinant Bovine Alpha-1-antiproteinase (SERPINA1)

<b>Product Code</b>	CSB-MP021053BO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P34955
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	GVLQGH AVQETDDTSH QEAACHKIAP NLANFAFSIY HHLAHQSNTS NIFFSPVIA SAFAMLSLGA KGNTHTTEILK GLGFNLTELA EAEIHKGFQH LLHTLNQPNH QLQLTTGNGL FINESAKLVD TFLEDVKNLY HSEAFSINFR DAEEAKKKIN DYVEKGSYGK IVELVKVLDP NTVFALVNYI SFKKGWKEKPF EMKHTTERDF HVDEQTTVKV PMMNRLGMFD LHYCDKLASW VLLLDYVGNV TACFILPDLG KLQQLLEDKLN NELLAKFLEK KYASSANLHL PKLSISETYD LKSVLGDVGI TEVFSRADL SGITKEQPLK VSKALHKAAL TIDEKGTEAV GSTFLEAIPM SLPPDVEFNR PFLCILYDRN TKSPLFVGKV VNPTQA
<b>Source</b>	Mammalian cell
<b>Target Names</b>	SERPINA1
<b>Protein Names</b>	Recommended name: Alpha-1-antiproteinase Alternative name(s): Alpha-1-antitrypsin Alpha-1-proteinase inhibitor Serpin A1
<b>Expression Region</b>	25-416
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This protein is secreted and is a serine protease inhibitor whose targets include elastase, plasmin, thrombin, trypsin, chymotrypsin, and plasminogen activator. Defects in this gene can cause emphysema or liver disease. Several transcript variants encoding the same protein have been found for this gene.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.