



Recombinant Bovine Alpha-1-antiproteinase (SERPINA1)

Product Code	CSB-YP021053BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P34955
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	≥85% (SDS-PAGE)
Sequence	GVLQGH AVQETDDTSH QEAACHKIAP NLANFAFSIY HHLAHQSNTS NIFFSPVSIA SAFAMLSLGA KGNTHTTEILK GLGFNLTELA EAEIHKGFQH LLHTLNQPNH QLQLTTGNGL FINESAKLVD TFLEDVKNLY HSEAFSINFR DAEEAKKKIN DYVEKGSBGK IVELVKVLDP NTVFALVNYI SFKKGKWEKPF EMKHTTERDF HVDEQTTVKV PMMNRLGMFD LHYCDKLASW VLLLDYVGNV TACFILPDLG KLQQLDKLN NELLAKFLEK KYASSANLHL PKLSISETYD LKSVLGDVGI TEVFSDRADL SGITKEQPLK VSKALHKAAL TIDEKGTEAV GSTFLEAIPM SLPPDVEFNR PFLCILYDRN TKSPLFVGKV VNPTQA
Source	Yeast
Target Names	SERPINA1
Protein Names	Recommended name: Alpha-1-antiproteinase Alternative name(s): Alpha-1-antitrypsin Alpha-1-proteinase inhibitor Serpin A1
Expression Region	25-416
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein is secreted and is a serine protease inhibitor whose targets include elastase, plasmin, thrombin, trypsin, chymotrypsin, and plasminogen activator. Defects in this gene can cause emphysema or liver disease. Several transcript variants encoding the same protein have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.