




Recombinant *Saccharomyces cerevisiae* TATA-box-binding protein (SPT15)

Product Code	CSB-YP023239SVG
Relevance	General transcription factor that functions at the core of the DNA-binding general transcription factor complex TFIID. Binding of TFIID to a promoter (with or without TATA element) is the initial step in preinitiation complex (PIC) formation. TFIID plays a key role in the regulation of gene expression by RNA polymerase II through different activities such as transcription activator interaction, core promoter recognition and selectivity, TFIIA and TFIIB interaction, chromatin modification (histone acetylation by TAF1), facilitation of DNA opening and initiation of transcription.
Abbreviation	Recombinant <i>Saccharomyces cerevisiae</i> SPT15 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P13393
Alias	TATA sequence-binding protein ;TBPTATA-binding factorTATA-box factorTranscription factor DTranscription initiation factor TFIID TBP subunit
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	Greater than ≥ 90% as determined by SDS-PAGE.
Sequence	ADEERLKEFKEANKIVFDPNTRQVWENQNRDGTKPATTFQSEEDIKRAAPESE KDTSATSGIVPTLQNIIVATVTLGCRLDLKTVLHARNAEYNPKRFAAVIMRIREP KTTALIFASGKMVVTGAKSEDDSKLASRKYARIIQKIGFAAKFTDFKIQNIVGSC DVKFPIRLEGLAFSHGTFSSYEPELFPGLIYRMVVKPKIVLLIFVSGKIVLTGAKQR EEIYQAFEAIYPVLSEFRKM
Research Area	Others
Source	Yeast
Target Names	SPT15
Protein Names	Recommended name: TATA-box-binding protein Alternative name(s): TATA sequence-binding protein Short name= TBP TATA-binding factor TATA-box factor Transcription factor D Transcription initiation factor TFIID TBP subunit
Expression Region	2-240aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	28.9kDa

**Protein Length**

Full Length of Mature Protein

Image

(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.