



Recombinant Phytolacca dioica Ribosome-inactivating protein PD-L1/PD-L2

Product Code	CSB-MP306294PJJ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P84853
Product Type	Recombinant Protein
Immunogen Species	Phytolacca dioica (Bella sombra tree) (Phytolacca arborea)
Purity	≥85% (SDS-PAGE)
Sequence	INTITYDAGN TTINKYATFM ESLRNEAKDP SLQCYGIPML PNNSSTIKYL LVKLQGASQK TITLMLRRNN LYVMGYSDPF NGNCRYHIFN DITGTERTNV ENTLCSSSSS RDAKPINYNS LYSTLEKKA E VNSRSQVQLG IQILSSDIGK ISGQSSFTDK TEAKFLLVAI QMVSEAAARFK YIENQVKTNF NRDFSPNDKI LDLEENWGKI STAIHDATNG ALPKPLELKN ADGTKWIVLR VDEIKPDMGL LNYVNGTCQT T
Source	Mammalian cell
Protein Names	Recommended name: Ribosome-inactivating protein PD-L1/PD-L2 EC=3.2.2.22 Alternative name(s): rRNA N-glycosidase PD-L1/PD-L2
Expression Region	1-261
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.