



Recombinant Mouse Pyruvate dehydrogenase E1 component subunit alpha, somatic form, mitochondrial (Pdha1)

Product Code	CSB-BP017715MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P35486
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	F ANDATFEIKK CDLHRLEEGP PVTTVLTRED GLKYYRMMQT VRRMELKADQ LYKQKIIRGF CHLCDGQEAC CVGLEAGINP TDHLITAYRA HGFTFTRGLP VRAILAEALTG RRGGCAKGGK GSMHMYAKNF YGGNGIVGAQ VPLGAGIALA CKYNGKDEVC LTLYGDGAAN QGQIFEAYNM AALWKLPCIF ICENNRGGMG TSVERRAAST DYYKRGDFIP GLRVDGMDIL CVREATKFAA AYCRSGKGPI LMELQTYRYH GHSMSDPGVS YRTREEIQEV RSKSDPIMLL KDRMVNSNLA SVEELKEIDV EVRKEIEDAA QFATADPEPP LEELGYHIYS SDPPFEVRGA NQWIKFKSVS
Source	Baculovirus
Target Names	Pdha1
Protein Names	Recommended name: Pyruvate dehydrogenase E1 component subunit alpha, somatic form, mitochondrial EC= 1.2.4.1 Alternative name(s): PDHE1-A type I
Expression Region	30-390
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.