




Recombinant Mycoplasma hyopneumoniae 46 kDa surface antigen (p46)

Product Code	CSB-YP314621MLQ
Abbreviation	Recombinant Mycoplasma hyopneumoniae 46KDA surface antigen protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0C0J8
Alias	p46
Product Type	Recombinant Protein
Immunogen Species	Mycoplasma hyopneumoniae (strain J / ATCC 25934 / NCTC 10110)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	CGQTESGSTSDSKPQAETLKHKVSNSIRIALTDPDNPRWISAQKDIISYVDET EAATSTITKNQDAQNNWLTQQANLSPAPKGFIIAPENSGSGVGTAVNTIADKGIPI VAYDRLITGSDKYDWYVSFDNEKVGELQGLSLAAGLLGKEDGAFDSIDQMNE YLKSHMPQETISFYTIAGSQDDNNSQYFYNGAMKVLKELMKNSQNKIIDLSPEG ENAVYVPGWNYGTAGQRIQSFLTINKDPAGGNKIKAVGSKPASIFKGFLAPND GMAEQAITKLKLEGFDTQKIFVTGQDYNDKAKTFIKDGDQNMITYKPKVVLGKV AVEVLRVLIAKKNKASRSEVENELKAKLPNISFKYDNQTYKVQGKNINTILVSPVI VTKANVDNPDA
Research Area	Others
Source	Yeast
Target Names	p46
Protein Names	Recommended name: 46 kDa surface antigen Alternative name(s): p46
Expression Region	28-416aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	44.5kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.