



# Recombinant Rat Thyroxine-binding globulin (Serpina7)

<b>Product Code</b>	CSB-MP021063RA
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P35577
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	APHNSSEGKV TTCHLPQQNA TLYKMPSINA DFAFRLYRKL SVENPDLNIF FSPVSISAAL AMLSFGSGSS TQTQILEVLG FNLTDTPVKE LQQGFQHLIC SLNFPNNELE LQMGNVAVFIG QQLKPLAKFL DDVKTLYETE VFSTDFSNVS AAQHEINSYV EKQTKGKIVG LIQDLKLNII MILVNYIHFQ AQWANPFRVS KTESSNFSV DKSTTVQVPM MHQLEQYYHY VDVELNCTVL QMDYSANALA LFVLPKEGHM EWVEAAMSSK TLKKWNHLLQ KGWVELFVPK FSISATYDLG STLQKMGMRD AFAESADFPG ITKDNGLKLS YAFHKAVLHI GEEGTKEGAS PEAGSLDQPE VAPLHAVIRL DRTFLLMILE KRTRSVLFLG KVVDPSTKE
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Serpina7
<b>Protein Names</b>	Recommended name: Thyroxine-binding globulin Alternative name(s): Serpin A7 T4-binding globulin
<b>Expression Region</b>	21-418
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.