



Recombinant Lipase

Product Code	CSB-MP321696RDA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P19515
Product Type	Recombinant Protein
Immunogen Species	Rhizomucor miehei
Purity	>85% (SDS-PAGE)
Sequence	SIDGGI RAATSQEINE LTYYTTLNAN SYCRTVIPGA TWDCIHCDAT EDLKIKTWS TLIYDTNAMV ARGDSEKTIY IVFRGSSSIR NWIADLTFVP VSYPVSGTK VHKGFLDSYG EVQNELVATV LDQFKQYPSY KVAVTGHSLG GATALLCALD LYQREEGLSS SNLFLYTQQQ PRVGDPAFAN YVYSTGIPYR RTVNERDIVP HLPPAAFGL HAGEEYWITD NSPETVQVCT SDLETSDCSN SIVPFTSVLD HLSYFGINTG LCT
Source	Mammalian cell
Protein Names	Recommended name: Lipase EC= 3.1.1.3 Alternative name(s): Triacylglycerol lipase
Expression Region	95-363
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4? for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.