




ureA Antibody

Product Code	CSB-PA14347A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14916
Immunogen	Recombinant Helicobacter pylori Urease subunit alpha protein (1-238AA)
Raised In	Rabbit
Species Reactivity	Helicobacter pylori
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Urease subunit alpha (EC 3.5.1.5) (Urea amidohydrolase subunit alpha), ureA, hpuA
Immunogen Species	Helicobacter pylori
Research Area	Others
Target Names	ureA

Image



Western blot
 All lanes: ureA antibody at 2µg/ml + Helicobacter pylori bacteria liquid
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 27 kDa
 Observed band size: 27 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.