




ACPP Antibody

Product Code	CSB-PA11309A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15309
Immunogen	Recombinant Human Prostatic acid phosphatase protein (38-386AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prostatic acid phosphatase (PAP) (EC 3.1.3.2) (5'-nucleotidase) (5'-NT) (EC 3.1.3.5) (Ecto-5'-nucleotidase) (Thiamine monophosphatase) (TMPase) [Cleaved into: PAPf39], ACPP
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	ACPP

Image



Immunohistochemistry of paraffin-embedded human prostate tissue using CSB-PA11309A0Rb at dilution of 1:20

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.