



Recombinant Human Lactadherin (MFGE8)

Product Code	CSB-MP013752HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q08431
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	LDICSKN PCHNGGLCEE ISQEVRGDVF PSYTCTCLKG YAGNHCETKC VEPLGMENGN IANSQIAASS VRVTFGLGLQH WPELARLNR AGMVNAWTPS SNDDNPWIQV NLLRRMWVTG VVTQGASRLA SHEYLKAFKV AYSLNGHEFD FIHDVNKKHK EFGVNWKNNA VHVNLFETPV EAQYVRLYPT SHTACTLRF ELLGCELNGC ANPLGLKNNS IPDKQITASS SYKTWGLHLF SWNPSYARLD KQGNFNAWVA GSYGNDQWLQ VDLGSSKEVT GIITQGARNF GSVQFVASYK VAYSNDSANW TEYQDPRTGS SKIFPGNWDN HSHKKNLFET PILARYVRIL PVAWHNRIL RLELLGC
Source	Mammalian cell
Target Names	MFGE8
Expression Region	24-387
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.