




ANXA3 Antibody

Product Code	CSB-PA001844YA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12429
Immunogen	Recombinant Human Annexin A3 protein (5-322AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:10-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Annexin A3 (35-alpha calcimedin) (Annexin III) (Annexin-3) (Inositol 1,2-cyclic phosphate 2-phosphohydrolase) (Lipocortin III) (Placental anticoagulant protein III) (PAP-III), ANXA3, ANX3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ANXA3

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA001844YA01HU at dilution of 1:20