




IL18 Antibody

Product Code	CSB-PA614514YA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14116
Immunogen	Recombinant Human Interleukin-18 protein (1-193AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-18 (IL-18) (Ibctadekin) (Interferon gamma-inducing factor) (IFN-gamma-inducing factor) (Interleukin-1 gamma) (IL-1 gamma), IL18, IGIF IL1F4
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL18

Image



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA614514YA01HU at dilution of 1:20

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.