




Col1a1 Antibody

Product Code	CSB-PA005727YA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02454
Immunogen	Recombinant Rat Collagen alpha-1(I) chain protein (955-1207AA)
Raised In	Rabbit
Species Reactivity	Rat, Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:10-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Others
Target Names	Col1a1

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA005727YA01RA at dilution of 1:20