




# Col1a1 Antibody

<b>Product Code</b>	CSB-PA005727YA01RA
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P02454
<b>Immunogen</b>	Recombinant Rat Collagen alpha-1(I) chain protein (955-1207AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Rat, Human
<b>Tested Applications</b>	ELISA, IHC; Recommended dilution: IHC:1:10-1:100
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Research Area</b>	Others
<b>Target Names</b>	Col1a1

## Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA005727YA01RA at dilution of 1:20

**Usage** For Research Use Only. Not for use in diagnostic or therapeutic procedures.