



IL18 Antibody, Biotin conjugated

Product Code	CSB-PA614514YD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14116
Immunogen	Recombinant Human Interleukin-18 protein (1-193AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-18 (IL-18) (Iboctadekin) (Interferon gamma-inducing factor) (IFN-gamma-inducing factor) (Interleukin-1 gamma) (IL-1 gamma), IL18, IGIF IL1F4
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL18
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.