



Col1a1 Antibody, FITC conjugated

Product Code	CSB-PA005727YC01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02454
Immunogen	Recombinant Rat Collagen alpha-1(I) chain protein (955-1207AA)
Raised In	Rabbit
Species Reactivity	Rat
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Collagen alpha-1(I) chain (Alpha-1 type I collagen), Col1a1
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Others
Target Names	Col1a1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.