



IktA Antibody, HRP conjugated

Product Code	CSB-PA314743EB01ESE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C085
Immunogen	Recombinant Mannheimia haemolytica Leukotoxin protein (715-953AA)
Raised In	Rabbit
Species Reactivity	Mannheimia haemolytica
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leukotoxin (Lkt), IktA
Immunogen Species	Mannheimia haemolytica
Research Area	Others
Target Names	IktA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.