



IktA Antibody

Product Code	CSB-PA314743EA01ESE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C085
Immunogen	Recombinant Mannheimia haemolytica Leukotoxin protein (715-953AA)
Raised In	Rabbit
Species Reactivity	Mannheimia haemolytica
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leukotoxin (Lkt), IktA
Immunogen Species	Mannheimia haemolytica
Research Area	Others
Target Names	IktA