





ANGPT1 Antibody

Product Code	CSB-PA05739A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15389
Immunogen	Recombinant Human Angiopoietin-1 protein (16-498AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Angiopoietin-1 (ANG-1), ANGPT1, KIAA0003
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	ANGPT1

Image



Western blot
All lanes: ANGPT1 antibody at 2µg/ml + A549 whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted bands size: 58 kDa
Observed bands size: 58 kDa



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA05739A0Rb at dilution of 1:50



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.