




Cd46 Antibody

Product Code	CSB-PA004939LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O88174
Immunogen	Recombinant Mouse Membrane cofactor protein (45-329AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Membrane cofactor protein (CD antigen CD46), Cd46, Mcp
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Cd46

Image



Western Blot

Positive WB detected in: Mouse skeletal muscle tissue


All lanes: Cd46 antibody at 4µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 41, 34 kDa

Observed band size: 41 kDa



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA004939LA01MO at dilution of 1:100