




glsA1 Antibody

Product Code	CSB-PA302868EA01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P77454
Immunogen	Recombinant Escherichia coli Glutaminase 1 protein (1-310AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutaminase 1 (EC 3.5.1.2), glsA1, ybaS
Immunogen Species	Escherichia coli
Research Area	Signal Transduction
Target Names	glsA1

Image



Western blot
 All lanes: glsA1 antibody at 2µg/ml + DH5a whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 33 kDa
 Observed band size: 33 kDa

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.