




Thbs3 Antibody

Product Code	CSB-PA023489EA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05895
Immunogen	Recombinant Mouse Thrombospondin-3 protein (24-350AA)
Raised In	Rabbit
Species Reactivity	Mouse, Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thrombospondin-3, Thbs3, Tsp3
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Thbs3

Image



Western blot
 All lanes: Thbs3 antibody at 16µg/ml + HepG2 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 105 kDa
 Observed band size: 105 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.