



Tlr7 Antibody, FITC conjugated

Product Code	CSB-PA023606EC01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P58681
Immunogen	Recombinant Mouse Toll-like receptor 7 protein (27-348AA)
Raised In	Rabbit
Species Reactivity	Mouse
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Toll-like receptor 7, Tlr7
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Tlr7