




Tlr7 Antibody

Product Code	CSB-PA023606EA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P58681
Immunogen	Recombinant Mouse Toll-like receptor 7 protein (27-348AA)
Raised In	Rabbit
Species Reactivity	Mouse, Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Toll-like receptor 7, Tlr7
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Tlr7

Image



Western Blot

Positive WB detected in: Jurkat whole cell lysate, K562 whole cell lysate, A549 whole cell lysate

All lanes: Tlr7 antibody at 2µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 122 kDa

Observed band size: 122 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.