



IMUP Antibody

Product Code	CSB-PA003450DA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9GZP8
Immunogen	Recombinant Human Immortalization up-regulated protein (1-85AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Immortalization up-regulated protein (Hepatocyte growth factor activator inhibitor type 2-related small protein) (H2RSP) (HAI-2-related small protein), IMUP, C19orf33
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	IMUP