



TNF Antibody, HRP conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA023955DB01SH |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P23383 |
| Immunogen | Recombinant Sheep Tumor necrosis factor protein (78-234AA) |
| Raised In | Rabbit |
| Species Reactivity | Sheep |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Tumor necrosis factor (Cachectin) (TNF-alpha) (Tumor necrosis factor ligand superfamily member 2) (TNF-a) [Cleaved into: Tumor necrosis factor, membrane form (N-terminal fragment) (NTF); Intracellular domain 1 (ICD1); Intracellular domain 2 (ICD2); C-domain 1; C-domain 2; Tumor necrosis factor, soluble form], TNF, TNFA TNFSF2 |
| Immunogen Species | Ovis aries (Sheep) |
| Research Area | Others |
| Target Names | TNF |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |