



PYC1 Antibody

Product Code	CSB-PA319972LA01SVG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11154
Immunogen	Recombinant Baker's yeast Pyruvate carboxylase 1 protein (18-470AA)
Raised In	Rabbit
Species Reactivity	Baker's yeast
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pyruvate carboxylase 1 (EC 6.4.1.1) (Pyruvic carboxylase 1) (PCB 1), PYC1, PYV
Immunogen Species	Baker's yeast
Research Area	Others
Target Names	PYC1