



IL1B Antibody, Biotin conjugated

Product Code	CSB-PA011614HD01SH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21621
Immunogen	Recombinant Sheep Interleukin-1 beta protein (114-266AA)
Raised In	Rabbit
Species Reactivity	Sheep
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-1 beta (IL-1 beta), IL1B
Immunogen Species	Ovis aries (Sheep)
Research Area	Others
Target Names	IL1B