




IL1B Antibody

Product Code	CSB-PA011614HA01SH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21621
Immunogen	Recombinant Sheep Interleukin-1 beta protein (114-266AA)
Raised In	Rabbit
Species Reactivity	Sheep
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-1 beta (IL-1 beta), IL1B
Immunogen Species	Ovis aries (Sheep)
Research Area	Others
Target Names	IL1B

Image



Western blot

All lanes: IL1B antibody at 2µg/ml

Lane 1: Recombinant Interleukin-1 beta protein 0.05µg

Lane 2: Recombinant Interleukin-1 beta protein 0.5µg

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 31 kDa

Observed band size: 31 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.