



KCNJ1 Antibody, HRP conjugated

Product Code	CSB-PA012047LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48048
Immunogen	Recombinant Human ATP-sensitive inward rectifier potassium channel 1 protein (178-391AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP-sensitive inward rectifier potassium channel 1 (ATP-regulated potassium channel ROM-K) (Inward rectifier K(+) channel Kir1.1) (Potassium channel, inwardly rectifying subfamily J member 1), KCNJ1, ROMK1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KCNJ1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.