




C5 Antibody

Product Code	CSB-PA003995YA01PI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01032
Immunogen	Recombinant Pig Complement C5a anaphylatoxin protein (1-74AA)
Raised In	Rabbit
Species Reactivity	Pig
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:10-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Complement C5a anaphylatoxin, C5
Immunogen Species	Pig
Research Area	Immunology
Target Names	C5

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA003995YA01PI at dilution of 1:10

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.