



C5 Antibody, HRP conjugated

Product Code	CSB-PA003995YB01PI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01032
Immunogen	Recombinant Pig Complement C5a anaphylatoxin protein (1-74AA)
Raised In	Rabbit
Species Reactivity	Pig
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Complement C5a anaphylatoxin, C5
Immunogen Species	Pig
Research Area	Immunology
Target Names	C5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.