



C5 Antibody, HRP conjugated

Product Code	CSB-PA003995YB01PI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01032
Immunogen	Recombinant Pig Complement C5a anaphylatoxin protein (1-74AA)
Raised In	Rabbit
Species Reactivity	Pig
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Complement C5a anaphylatoxin, C5
Immunogen Species	Pig
Research Area	Immunology
Target Names	C5