



Sod1 Antibody, Biotin conjugated

Product Code	CSB-PA022397LD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08228
Immunogen	Recombinant Mouse Superoxide dismutase [Cu-Zn] protein (2-154AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Superoxide dismutase [Cu-Zn] (EC 1.15.1.1), Sod1
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Sod1