



glsA1 Antibody, FITC conjugated

Product Code	CSB-PA302868EC01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P77454
Immunogen	Recombinant Escherichia coli Glutaminase 1 protein (1-310AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutaminase 1 (EC 3.5.1.2), glsA1, ybaS
Immunogen Species	Escherichia coli
Research Area	Signal Transduction
Target Names	glsA1